



DermaRed

LED Age-Defying Therapy

Reverse Signs of Prematurely Aging Skin

Sound bite: "Shine the light on firmer, younger-looking skin."

What is it?

4-in-1 youth-activating system provides four technologically superior treatments into one tool to minimize the appearance of lines and wrinkles, improve skin tone and elasticity, deep cleanse pores, and infuse water-based treatments for firmer, smoother, more lifted and radiant-looking skin. LED **red light**, ultrasonic waves, and positive and negative ions visibly help reverse signs of prematurely aging skin after every treatment.



What does it do?

- **LED red light** airbrushes away the appearance of lines & wrinkles; improves the appearance of age spots & photo damage
- **Ultrasonic waves** stimulate blood circulation & improve the appearance of skin tone & texture
- **Positive ions** deep cleanse pores of excess oil & impurities
- **Negative ions** accelerate the absorption of serums & other water-based skincare products

How do I use it?

- Apply water-based product, serum, or gel to targeted areas.
- Press power button on, then ION+ button, and move unit across face to cleanse pores.
- Press LED button, SONIC button, and ION- button, and move unit across targeted areas in circular motions, keeping skin moist. Do not use near the eyes or eyelids.

What skin types and conditions does it benefit?

Skin types: All skin types

Skin conditions: Lines & wrinkles, uneven tone, enlarged pores, photo damage

Caution: Keep out of reach of children. Operate this unit only according to its intended use as described in the product manual. Do not use this unit if you have serious allergies, infected skin, or a heart condition. Do not immerse the unit in water. Keep the unit free of lint and hair.