

## TOOL

# BioSculpt7 LED Mask

Neck & Décolleté

**sound bite:** "7-in-1 skin-rejuvenating targeted light therapy"

FDA 510(k)-cleared multi-color LED neck mask gently smooths, tones, and revitalizes the appearance of crepey, lined, or flaccid-looking skin on the neck and décolleté. Featuring 7 LED light options for customized sessions, the flexible design makes advanced neck and décolleté treatments effortless and comfortable.

### benefits

Improves appearance of sagging neck and décolleté

Targets the look of crepey, lined skin

Encourages a visibly lifted appearance

Brightens the appearance of dull, tired skin

### light modes, wavelengths (nm) & targeted benefits

- Red Light** (630-640nm) Helps improve appearance of firmness & smoothness
- Green Light** (520-530nm) Supports the appearance of uneven skin tone
- Blue Light** (470-480nm) Helps promote clearer-looking skin
- Cyan Light** (520-530nm) Helps soothe irritated skin
- Yellow Light** (590-600nm) Brightens appearance of dull, tired skin
- Purple Light** (380-390nm) Revitalizes appearance of fatigued, lackluster skin
- White Light** (780-790nm) Supports overall skin brightness & vitality

### directions

Connect fully-charged controller to mask cable. Apply mask on CLEAN, DRY SKIN with LED lights facing skin. Use 10-15 minutes daily for best results. For full directions, scan QR code.

Skin may appear red immediately after treatment, usually returning to normal within 1-2 hours. Some users may notice pigmentation changes, flaking, peeling, itchiness, or dryness, which typically resolve within a few days. Temporary headaches may occur and can be relieved with rehydration or an over-the-counter analgesic.

**! contraindications** For external use only. Avoid direct eye exposure. Do not use over broken, irritated, or inflamed skin. Avoid using directly over the thyroid in cases of thyroid disease. Individuals with photosensitive epilepsy or taking photosensitizing medications should consult a physician before use. Discontinue use if redness, discomfort, or irritation occurs. Keep out of reach of children. Consult a physician before use if you are pregnant or have a serious medical condition. Do not use if mask, cable, or controller is damaged.



### skin type

- ☒ Oily
- ☒ Normal to Oily
- ☒ Normal
- ☒ Normal to Dry
- ☒ Dry

### skin concern

- ☒ Fine lines
- ☒ Deep lines
- ☒ Loss of elasticity
- ☒ Hyperpigmentation
- ☒ Blemishes