

Oxygenated Mist for Face and Body



SOUND BITE: "Help skin breathe, from head to toe."

Problem: Every day, oxygen levels in skin decline, bombarded by environmental stressors that trigger dehydration, wrinkles, sagging, uneven skin tone, blemishes, and inflammation.

Solution: **O₂ Infuser** dispenses oxygenated skincare nutrients in an ultra-fine mist that instantly hydrates, nourishes, and brightens skin on face and body. Oxygenated, micronized mist propels skincare product molecules into skin for an immediate dewiness and healthy glow. Easy-to-use dispenser leaves skin cool, refreshed, and rejuvenated.

FOR USE WITH WATER-BASED PRODUCTS ONLY.

WHAT DOES IT DO?

- Saturates skin in oxygenated skincare mist
- Instantly hydrates and creates a dewy glow
- Cools, refreshes, and rejuvenates face AND body
- Helps ward off effects of environmental stressors

DIRECTIONS

Pour Oh My Glow Activator into mini-flask attached to top of O₂ Infuser. Add Oh My Glow Serum. Close lid and gently shake flask to blend liquid. Press power button once for high intensity or twice for low intensity. Pull back on sprayer switch to dispense oxygenated mist. Press power button three times to turn off machine.

SKIN TYPES

all skin types

CONDITIONS

Dryness, lines & wrinkles, uneven tone

CAUTION

For use with water-based products only. Do not use when unit is plugged in. Do not use while unit is charging. Do not leave product inside unit for a prolonged period of time. Operate only according to intended use as described in user manual. Keep out of reach of children.

5/11/21