

Oxygenated Mist for Face and Body



SOUND BITE: "Help skin breathe, from head to toe."

Problem: Every day, oxygen levels in skin decline, bombarded by environmental stressors that trigger dehydration, wrinkles, sagging, uneven skin tone, blemishes, and inflammation.

Solution: O₂ INFUSER dispenses oxygenated skincare nutrients in an ultra-fine mist that instantly hydrates, nourishes, and brightens skin on face and body. Oxygenated, micronized mist propels skincare product molecules into skin for an immediate dewiness and healthy glow. Easy-to-use dispenser leaves skin cool, refreshed, and rejuvenated.

For use with water-based products only.

WHAT DOES IT DO?

- Saturates skin in oxygenated skincare mist
- Instantly hydrates & creates a dewy glow
- Cools, refreshes & rejuvenates face AND body
- Helps ward off effects of environmental stressors

DIRECTIONS

Attach flask to infuser. Fill halfway with **Oh My Glow Activator**. Add 2-3 drops **Oh My Glow Serum** (or more for concentrated treatment). Close lid. Shake flask to blend. Press power button once for high intensity or twice for low intensity. Pull sprayer switch to dispense oxygenated mist.

SKIN TYPES

all skin types

CONDITIONS

Dryness, lines & wrinkles, uneven tone

CAUTION

For use with water-based products only. Do not use when unit is plugged in. Do not use while unit is charging. Do not leave product inside unit for a prolonged period of time. Operate only according to intended use as described in user manual. Keep out of reach of children.