

TOOL

Gua Sha ThermoLift

Lift & Contour

sound bite: “Next-gen sculpting tool.”

Gua sha-inspired, multi-functional sculpting tool delivers micro vibrations, red and blue LED light, and therapeutic heat to boost microcirculation, reduce puffiness, and soften lines for a smoother, more lifted look.

Red LED targets signs of aging, blue LED addresses breakouts and visible inflammation—especially in skin exposed to daily environmental stressors—and thermal warmth enhances product absorption for maximum results.

what does it do?

LIFT: Microvibrations & gua sha strokes help tighten & firm

CONTOUR: Microvibrations & gua sha strokes help sculpt & tone

DETOX: Blue light addresses breakouts, impurities & inflammation

INFUSE: Red light & thermal warmth enhance product absorption

directions

After charging tool, apply Le Mieux skincare product. Press and hold Power On button. Select mode, then glide the tool gently in 1 direction along Meridian lines for 5-10 minutes using Meridian points as reference. Ensure the heating function feels comfortable on skin. Discontinue use if uncomfortable. For best results, advanced education on gua sha massage is recommended.

contraindications

This unit should not be used by pregnant women or individuals with a pacemaker or other implanted electronic medical devices. If client experiences any discomfort, irritation, or adverse reactions, discontinue use immediately and consult a healthcare professional. Operate only according to intended use as described in user guide. Keep out of reach of children.



skin type

- ☒ Oily
- ☒ Normal to Oily
- ☒ Normal
- ☒ Normal to Dry
- ☒ Dry

skin concerns

- ☒ Fine lines
- ☒ Loss of elasticity
- ☒ Sensitivity
- ☒ Blemishes