

TOOL

Ionized Oxygen Infuser

Mist & Glow

sound bite: *"Instant Line-Filler."*

Easy-to-use, customizable infuser dispenses an ultra-fine mist to instantly hydrate, plump fine lines, volumize, and create a dewy glow.

Use with Oh My Glow Activator and Oh My Glow Serum, or customize with skin-specific Le Mieux toners and water-based serums.

what does it do?

Saturates skin with ultra-fine mist to fill lines & volumize

Instantly hydrates & creates a dewy glow

Cools, refreshes & rejuvenates face & body

Customize with skin-specific Le Mieux toners and water-based serums

directions

Pour Oh My Glow Activator into mini-flask attached to top of infuser. Add Oh My Glow Serum. Close lid and gently shake flask to blend liquid. Press power button and apply oxygenated mist. Press power button to turn off machine.

contraindications

For use with water-based products only. Do not use when unit is plugged in. Do not use while unit is charging. Do not leave product inside unit for a prolonged period of time. Operate only according to intended use as described in user manual. Keep out of reach of children.



skin type

- ☒ Oily
- ☒ Balanced to Oily
- ☒ Balanced
- ☒ Balanced to Dry
- ☒ Dry

skin concern

- ☒ Fine lines
- ☒ Dehydration
- ☒ Sensitivity